



***Brucella* spp. Specimen Selection**

+ Serum

- The diagnosis of brucellosis is frequently achieved by serology. An acute & convalescent phase specimen should be collected (21 days apart)

+ Blood or bone marrow

- Brucellae most often isolated from these

+ Tissue (spleen, liver)

- Brucellae occasionally isolated



***Brucella* spp. Biosafety Alert**

- + Brucellosis is THE most commonly reported laboratory-associated bacterial infection.**
- + Cases have occurred in clinical laboratory settings by “sniffing” cultures, direct skin contact with cultures, and aerosol generating procedures**



Level A Laboratory Tests ***Brucella spp.***

- + Colonial morphology on SBA**
- + Gram stain morphology**
- + Oxidase**
- + Urea hydrolysis**
- + BSL 3**



***Brucella* spp. Key Level A Lab Tests**

- + Colonial morphology on SBA**
 - Fastidious
 - Visible growth may take 48 - 72 hrs
 - Small (0.5-1.0mm), convex, glistening
 - Non-hemolytic and non-pigmented



***Brucella* spp. Key Level A Lab Tests**

+ Gram Stain Morphology


- Tiny (very) - 0.5 - 0.7? x 0.6 - 1.5?
- Faintly staining
- Gram-negative coccobacilli

+ Technical Hints

- *B. melitensis*: the most coccal (rarely >1?)
- *B. abortus*: may be rod-like (2 - 3 ?)



***B. abortus*, gram stain, (x2400)**





B. abortus, gram stain, (x3200)





Brucella spp. **Key Level A Lab Tests**

+ Oxidase-positive

- *B. melitensis* (100%)
- *B. abortus* (96%)
- *B. suis* (95%)
- *B. canis* (72%)

+ Urea hydrolysis-positive

- *B. suis* & *B. canis* ~15 min
- *B. abortus* & *B. melitensis* ~24hr



Brucella spp. **Other Ox+/Ur+ GN Oxidizers**

- + **Achromobacter** grp B
- + **Acidovorax** spp
- + **Agrobacterium** spp
- + **EO-2/EO-3**
- + **Flavobacterium** spp
- + **Methylobacterium** spp
- + **Ochrobactrum**
- + **Pseudomonas** spp
- + **Riemerella**
- + **Roseomonas** spp
- + **O-2**
- + **II-i**



Brucella spp. **Rapid Method Results**

+ Dade Behring

- **MicroScan Conventional Panels:**
 - Currently, not in data base
- **MicroScan Rapid Panels:**
 - Currently, not in data base

+ Biolog

- Currently, not in data base



Brucella spp. **Rapid Method Results**

+ bioMérieux

- API 20E: Currently, not in data base.
 - At 48 hrs with only oxidase +, “Low discrimination” with a note: possible *Brucella* spp.



Brucella spp. **Rapid Method Results**

Vitek GNI + Card: Not in data base.

- ***B. abortus***: Nonfermenting GNB (asaccharolytic) or *Chryseobacterium indologenes*/*Brevundimonas vesicularis* or Questionable Biopattern
- ***B. canis***: *Acinetobacter* urea or Nonfermenting GNB/*Acinetobacter* urea



Brucella spp. **Rapid Method Results**

+ bioMérieux (cont'd)

- **Vitek GNI + Card: Not in data base**
 - *B. melitensis*: Nonfermenting GNB (asaccharolytic)
 - *Flavobacterium odoratum* GCMS (Good Confidence Marginal Separation)
 - *B. suis*: *Ochrobactrum anthropi* or unidentified
- **Vitek 2 GNB Card: IN THE DATA BASE**
 - *Brucella spp.*: True accuracy not yet determined



Brucella spp. **Technical Hints**

- + Misidentified as *Moraxella sp.***
 - Clin Inf Dis 1993; 17:1068-9
- + Reported as gram-positive cocci, mistaken for “slow-growing” *Staphylococcus sp.***



Brucella spp. **Review of Key Tests**

- + Tiny, faintly staining, gram-negative coccobacilli from blood or bone marrow
- + Slow growth on SBA, 2-3 days for colony appearance
- + Oxidase +
- + Urease +



Clostridium botulinum

Botulism



BOTULISM

- + The microbe, *C. botulinum*, was first described in 1897 by E. van Ermengem following a foodborne outbreak in Belgium
- + Botulism is a neuromuscular illness resulting from the action of a potent toxin



BOTULISM

- + The diagnosis of botulism is made clinically, i.e., based on the patient's case history and physical findings
- + Health care providers suspecting botulism should contact their State Health Department Epidemiology Office




BOTULISM

- + The Foodborne and Diarrheal Diseases Branch of CDC provides emergency consultation and support to public health authorities
- + Days: 404-639-2206, M-F, 0830-1630 EST
- + Nights/weekends: 404-639-2888



Laboratory Capacity for Botulinum Toxin Testing





BOTULISM

- + There are seven types of *botulinum* toxin:
 - Types A, B, E and in rare cases, F causes disease in humans
 - Types C and D cause disease in birds and mammals
 - Type G, identified in 1970, has yet to be confirmed as a cause of illness in man or animals



FOODBORNE BOTULISM

- + Infective dose: 0.001 µg/kg (type A, LD₅₀)
- + Incubation period: 18 - 36 hr (6hr to 10 d)
- + Dry mouth, double vision, droopy eyelids, dilated pupils
- + Generalized, progressive descending bilateral muscle weakness & paralysis
- + Respiratory failure and death
- + Mortality 5-10%, up to 25%



FOODBORNE BOTULISM

- + Among 309 persons with clinically diagnosed botulism reported to CDC from 1975 to 1988:**
 - Stool cultures for *C. botulinum*: 51% +**
 - Serum botulinum toxin testing: 37% +**
 - Stool botulinum toxin testing: 23% +**
- + Overall, at least one of the above tests was positive for 65% of all patients**



Level A Procedures for Botulism Event

- + Properly collected specimens are to be referred to designated testing laboratories**
- + Prior to the shipment of any botulism-associated specimen, the designated laboratory must be notified and approved by the State Health Epidemiology Office**



Level A Procedures for Botulism Event

- + Clinical specimens to be collected:**
 - 1. Serum**
 - 2. Gastric contents or vomitus**
 - 3. Feces or return from sterile water enema**
 - 4. Food samples**
- + Autopsy specimens:**
 - 1. Serum**
 - 2. Gastric and intestinal contents**



Botulism Biosafety Alert

- + Botulism toxins are extremely poisonous
- + Minute quantities acquired by ingestion, inhalation, or by absorption can cause profound intoxication and death
- + All materials suspected of containing toxin must be handled with CAUTION!



Botulism Biosafety Alert

- + Materials suspected of containing botulism toxin must be handled:
 - **Biological Safety Cabinet (Class II)**
 - **Laboratory Coats**
 - **Disposable surgical gloves**
 - **Face shield (as needed)**



Botulism Clinical Hints

+ The specimen of choice for confirming the diagnosis of botulism:

- 1. Serum toxicity testing**
- 2. Cultures* of debrided tissue from wound**
- 3. Cultures* of stool or incriminated food(s)**

***confirm isolates are toxigenic**



Bacillus anthracis

Anthrax



ANTHRAX

+ Three forms of human anthrax occur:

1. Cutaneous

2. Gastrointestinal

- Oropharyngeal
- Abdominal

3. Inhalation (Woolsorter's Disease)



Cutaneous anthrax

Vesicle development, day 2



Eschar formation, day 4





Anthrax Lesion on Neck





Inhalation Anthrax

- + Infective dose = 8,000 - 15,000 spores**
- + Incubation period = 1-6 days**
- + Duration of illness = 3-5 days**
- + Fever, malaise, and fatigue**
- + Short period of improvement = up to 2 days**
- + Abrupt respiratory distress...death <24hrs**
- + Person to person transmission = no**




Inhalation Anthrax

- + Chest x-ray is often pathognomonic: widened mediastinum with pleural effusions, without infiltrates.**
- + Late in the course of infection, gram stain of blood may reveal broad gram-positive rods.**
- + Only vegetative, encapsulated bacilli are present in tissues during infection, typically spores are not seen**



Mediastinal LN, microcolonies of *B anthracis*, Giemsa stain





Laboratory Handling of *Bacillus anthracis*, BSL 2

- + Anthrax vaccination is NOT required**
- + Laboratory coats and gloves must be worn**
- + Safety glasses/shields are recommended**
- + No hand/finger contact with mucosal surfaces**
- + Hand washing before leaving lab area**



Anthrax: Specimen Selection

- + Inhalation: Sputum and Blood**
- + Cutaneous: Vesicles and Eschar**
- + Gastrointestinal: Stool and Blood**



Bacillus anthracis **Key Level A Lab Tests**

- + Gram stain**
- + Growth characteristics on agar**
- + Sporulation, in air**
- + Motility**
- + (Penicillin-zone)***
- + (Capsule by India Ink)***
 - *optional**



Bacillus anthracis **Gram Stain Morphology**

- + Broad gram-positive rod: 1-1.5 X 3-5 μ**
- + Oval, central - subterminal spores: 1 X 1.5 μ
with no significant swelling of cell**
- + Spores are NOT usually present in clinical
specimens unless exposed to atmospheric
O₂**



***B. anthracis*, gram stain demonstrating spores**



4




B. anthracis **Colonial Morphology**

- + Colonial morphology of 18-24hr @ 35 C:**
 - Well isolated colonies are 2-5 mm in diameter**
 - Flat or slightly convex, irregularly round**
 - Edges: slightly undulate, often curly tailing edges**
 - Ground glass appearance**
 - “Sticky” consistency....stands up like beaten egg whites**




***B. anthracis*, colony on SBA**





“STICKY” consistency of *B. anthracis*’ colony on SBA





***B. anthracis* on blood (left) and bicarbonate agar (right) plate
44 (demonstrating capsule formation)**



Bacillus anthracis **Presumptive Identification**

- + Gram-positive, broad rod, catalase-positive, spore-positive, aerobe:
Bacillus sp.
- + Spores are oval and nonswelling with ground glass colony appearance:
Bacillus morphology group 1, includes
B. anthracis, *B. cereus*, *B. cereus var mycoides*, and *B. thuringiensis*



Bacillus anthracis **Presumptive Identification, con't**

- + Nonmotile: *B anthracis* and *B cereus* var *mycoides* (and *B. megaterium*)**
- + Nonhemolytic, penicillin-zone (15-20 mm) and forms capsule: Presumptive *B. anthracis***




CDC Laboratory Pearl

- + The most common *Bacillus sp.* submitted to CDC to r/o *B. anthracis* are non-motile *B. megaterium*
- + A non-motile, *Bacillus sp.*, is recovered from a blood culture. What is one of the fastest procedures to r/o anthrax ?



Bacillus anthracis **TECHNICAL HINTS**


- + Direct smears from clinical specimens:** encapsulated broad rods in short chains, 2 - 4 cells. India Ink will demonstrate capsule (gram stain will not)
- + Smears from SBA:** non-encapsulated broad rods in long chains






Packaging and Labeling of Infectious Substances

For all shipping and referrals:





**Infectious
substance label
for shipping all
specimens
and cultures**



Level A Laboratory Packaging/Shipping

- + Any clinical specimen (swabs, scrapings, body fluids, tissue) and/or “suspected” organisms that are referred to a higher level laboratory (B,C, or D) to rule out any of the agents discussed are to be classified as “INFECTIOUS SUBSTANCES”, and packaged and shipped accordingly.



Level A Laboratory Packaging/Shipping

- + All packaging must be UN certified for infectious substances, designated as: UN 4G/CLASS 6.2
- + NOTE: The above packaging is a “combination packaging” and must be used as a unit.



Level A Laboratory Packaging/Shipping

+ Primary Receptacle: WATERTIGHT

- Includes: glass, plastic, metal or screw-cap tubes.
- Label: Specimen ID label
- To ensure a leak-proof seal, screw caps must be fastened with tape, shrink seals or other comparable packaging



Level A Laboratory Packaging/Shipping

- Absorbent material must surround each primary receptacle and be sufficient to absorb the contents of primary container



Level A Laboratory Packaging/Shipping

- + Secondary Packaging: WATERTIGHT**
 - Primary receptacle must be placed within the secondary packaging surrounded by absorbent material
 - Itemized list of contents, e.g., *E. coli* on 5ml slant, affixed to secondary packaging
 - Secondary package must be also be labeled with the following information: name, address, and telephone number of shipper



Level A Laboratory Packaging/Shipping

- + Outer Packaging: ADEQUATE STRENGTH**
 - Secondary packaging (along with proper forms required from higher level laboratories) must be placed within the outer packaging unit**



Level A Laboratory Packaging/Shipping

- **Labels required on outer packaging:**
 - **Address label: higher level lab's address and phone number**
 - **Shipper's name, address and telephone number**
 - **Infectious substances label**
 - **Completed UN 2814 label**



Level A Laboratory Packaging/Shipping

- + Overpack: ADEQUATE STRENGTH**
 - A must with dry ice or cold packs
 - **Dry ice must not be packed inside the combination packaging**
 - Outer packaging (along with proper forms required from higher level laboratories) must be placed within the overpack
 - Labels required for overpack (see next slide)



Level A Laboratory Packaging/Shipping

+ Overpack: ADEQUATE STRENGTH

– Labels/forms required on overpack include:

- » **Higher level lab's address and phone #**
- » **Shipper's name, address and telephone number**
- » **Infectious substances label**
- » **Completed UN 2814 label**
- » **"Inner packages comply with prescribed specifications"**
- » **Shipper's Declaration for Dangerous Goods form**
- » **Label for ice/dry ice (if used)**
- » **Contact shipping company for additional**



Level A Laboratory Packaging/Shipping

SHIPPING CARRIER CRITERIA:

- + 1. Overnight delivery
- + 2. Package tracking capability
 - Such as: Airway Bill #



Websites of Interest

- + **Biosafety in Microbiology Lab**
 - www.cdc.gov/od/ohs/
- + **Guideline for Isolation Precautions**
 - www.cdc.gov/ncidod/hip/
- + **Public Health Image Library**
 - phil.cdc.gov/
- + **WHO: zoonotic diseases**
 - www.who.int/emc/