

Current expectations for laboratory testing and adverse smallpox vaccine reactions

- **Poxvirus Section (DVRD/NCID/CDC)**
- **World Health Collaborating Center for Smallpox and other Poxvirus Infections**

Inger Damon: Chief
Russell Regnery

Draft 1/20/03 Do not reproduce



Poxvirus intro:

2 Subfamilies:

Chordopoxvirinae (vertebrate poxviruses)

- Orthopoxvirus (variola, vaccinia, cowpox, monkeypox, raccoonpox, camelpox, skunkpox, volepox, ectromelia, taterapox)
- Parapoxvirus (orf, pseudocowpox, ...)
- Avipoxvirus (canarypox, fowlpox...)
- Capripoxvirus (goatpox, lumpy skin disease...)
- Leporipoxvirus (myxoma, fibroma...)
- Molluscipoxvirus (molluscum contagiosum)
- Yatapoxvirus (tanapox, Yaba monkey tumor)
- Entomopoxviridae (insect poxviruses)


Orthopoxviruses (including vaccinia and variola)

- 180 to 200 kbp double stranded DNA genomes with over 180 predicted ORFs
 - Encodes transcription and replication enzymes
 - Encodes multiple proteins aimed at evasion of immune defense molecules:
 - Soluble cytokine and chemokine binding proteins
 - Caspase inhibitors, inhibitors of protein kinase PKR
 - Cytoplasmic replication
 - Infectious forms: IMV, CEV, EEV
 - No known unique viral receptor protein
 - Host ranges vary
 - Variola (host specific...man) vs vaccinia (wide host range)
 - 350 X 270 nm brick shaped particles by cryoelectron microscopy
 - Antigenically similar; serologic cross reactivity

Orthopoxviruses: spectrum of human disease in normal host

- **Vaccinia**, cowpox: localized infection
 - **Variola**, monkeypox : systemic illness
-
- Camelpox: no known infections of humans
 - Ectromelia: no infection of humans
 - Raccoonpox, volepox, skunkpox: rare if any infection of humans with these animal pathogens

Febrile, vesicular rash illness algorithm for evaluating patients for smallpox



Draft 1/20/03 Do not reproduce

Key concepts of algorithm for evaluating pts for smallpox

- Product of practical experience and partnerships
- Uses existing resources for dx and exclusion of smallpox look-a-likes (especially 1.5m cases of chickenpox)
- Provide significant diagnostic benefits, even in absence of smallpox...encourages careful dx of other rash illnesses
- Minimize number of cases that require intensive investigation: focus attention where it is justified
- Provide rapid, thorough response to highly suspect cases
- Respond with vaccination if/when diagnosis confirmed

Draft 1/20/03 Do not reproduce

Smallpox dx: the Bottom Line

- To rapidly recognize and respond to the first case of smallpox quickly *without*
 - Generating high numbers of false alarms
 - Disrupting the health care and public health systems
 - Increasing public anxiety
- **Secondarily**, benefit health care (in absence of smallpox) by enhancing rash illness dx capability

Draft 1/20/03 Do not reproduce

Smallpox dx rational (pre-event)

- *Multiple samples (& multiple tests) to improve predictive value for positive variola lab dx in absence of disease in nature.*
- Multiple samples (e.g., biopsy) may provide means for dx of smallpox look-alikes; improves Public Health understanding of what can be clinically confused with smallpox

Vaccinia identification: lab expectations/considerations (I)

- **Most AE cases linked directly to vaccination history or contact of vaccinee.**
- **Vaccinia therapeutic options limited to VIG or Cidofovir (?)**
- **Decision to use VIG or Cidofovir made at clinical level**

Vaccinia identification: lab expectations/considerations (II)

- Identification of vaccinia in AE's will help improve Public Health understanding of AE's and vaccination risks; hopefully lead to better vaccines / vaccination strategies. Not a STAT function.
- LRN labs have means to detect vaccinia
 - But currently R-T PCR test considered by the FDA an “investigational device” or a “presumptive screening assay”.
 - Test results for pt management must be confirmed at CDC.
- Rule out other possible etiologies

Draft 1/20/03 Do not reproduce



Smallpox vs vaccinia: Lab tests may be similar but expectations for results and responses different

- Pre-event *smallpox* dx implies international crisis & mass vaccination
- *Smallpox* dx implies pt isolation and vaccination of contacts
- *Vaccinia* AE's expected to occur in small numbers (not the prelude to a pandemic!)
- AE patient care decisions based primarily on pt history and clinical considerations

Smallpox & generalized vaccinia: Differential Diagnosis & other rash illnesses

- **Varicella****
- **Disseminated herpes zoster***
- **Impetigo**
- **Drug eruptions**
- **Contact dermatitis**
- **Erythema multiforme**
- **Enteroviral infections**
- **Disseminated herpes simplex**
- **Scabies, insect bites**
- **Molluscum contagiosum (immunocompromised)**

Laboratory Testing to Rule Out Other Rash Causing Diseases

- VZV: DFA, PCR, EM, Immunohistochemistry
- Streptococcus, staphylococcus: Gram stain, rapid tests
- Drug eruptions, allergic dermatitis: skin biopsy, pathology
- Enterovirus infections: PCR, immunohistochemistry

Laboratory Testing to Rule Out Other Rash Causing Diseases - Continued

- Herpes simplex: PCR, EM, Immunohistochemistry and Culture (with caution)
- Scabies: Evidence of organisms
- Others as indicated from clinical impression:
 - *Rickettsia*: PCR
 - Syphilis

Varicella Rule-Out

- History and Clinical Findings
- Laboratory:
 - VZV DFA (two commercial products: Meridian Diagnostics and Chemicon)
 - PCR (“Homebrew”) plus LRN supported test
 - Electron microscopy
 - Refer All specimens with questionable results to DPH Laboratory for Confirmation

Generalized vaccinia, varicella, and smallpox: typical clinical features

	Generalized Vaccinia	Varicella (chickenpox)	Smallpox
Prodrome	Rare ??	Rare	Always
Distribution	Vac site + disseminated	Centripetal	Centrifugal
Lesions	Vac site + deep- seated pustules	Superficial vesicles	Deep-seated pustules
Progression	??	Rapid: some crusts <24 h	Slow: each stage 1-2 d
Stages of Dev't	Same (following vac) ?	Different	Same
Palms/soles	??	Extremely rare	Typical
Vac status	Recent smallpox vac or contact	Rare with VZV vac	Very rare with smallpox vac
Toxic	No	No	Yes

Smallpox (plus vaccinia?) Lab Algorithm

Low & Moderate Risk Specimens
(Green & Yellow Box)

Level A (hospitals)* and/or
LRN B/C Laboratories

VZV, HSV DFA

VZV, HSV, Enterovirus PCR
(where available)

Other appropriate diagnostic
assays, including viral culture

Positive = Diagnosis
of non-variola rash
illness

Negative:

- Re-evaluate patient condition and assess need for dermatologic and histologic testing, including tests for erythema multiforme.
- Obtain detailed information about possible **exposure to smallpox vaccines**.

Tests for vaccinia and generic orthopox may be indicated. IF vaccinia -, generic orthopox +, contact CDC immediately.

- If patient symptoms progress to more closely resemble smallpox, all specimens should be referred to CDC (and regional labs when needed).



Orthopoxvirus-specific assays for identification of:

- **Vaccinia (e.g., generalized vaccinia)**
- **Variola (smallpox)**

Lab methods for confirmation of orthopoxvirus dx

- PCR related methods for DNA identification, e.g., real-time PCR
- Electron microscopy
- Histopathology
- Culture
- Serology
 - Antigen detection (IFA/DFA, EIA ag capture)...not yet available/licensed

Draft 1/20/03 Do not reproduce

Vaccinia-related disease: ***potential* samples and tests**

● **Samples:**

- Lesion 'roofs' and crusts
- Vesicular fluids:
 - touch prep
 - Electron microscope grid
- Biopsy, autopsy

● **Tests:**

- R-T PCR (LRN)
- Electron microscopy
- Isolation in cell culture
- Antigen Capture & DFA (in the future)
- Histopath

Sample requirements for Poxvirus *DNA* identification

- ✓ Lesion 'roofs' and crusts
- ✓ Vesicular fluids (touch prep)
- ✓ Biopsy, PM autopsy
- ✓ Others (e.g. CSF?)

Sample requirements for *EM* poxvirus identification**

- ✓ Lesion 'roofs' and crusts
(grids subsequently made at lab)
- ✓ Vesicular fluids:
 - touch prep slide (grids made at lab)
 - electron microscope grids made on-site

**** requires access to EM**

Draft 1/20/03 Do not reproduce

Sample requirements for *vaccinia* (not *variola**) *isolation*

- ✓ Lesion 'roofs' and crusts
- ✓ Vesicular fluids:
 - touch prep slide (reconstituted at lab)
- ✓ Frozen biopsy including PM tissue (?)
- ✓ **BSL-II conditions; vaccinated lab workers preferable**

*Variola referred to CDC (pre-event)

Draft 1/20/03 Do not reproduce



Sample requirements for *histopathology*

✓ **Biopsy or autopsy
formalin fixed (not frozen)**

**(remember to save fresh frozen
bisected or duplicate sample for
isolation)**



Specimen collection

- Vaccinia and variola specimen collection essentially the same.
- *Check CDC website for recent updates in orthopox specimen collections specifics:*
 - <http://www.bt.cdc.gov/agent/smallpox/response-plan/files/guide-d.pdf>
 - Expect vaccinia specific updates in near future (first week December 2002)

Draft 1/20/03 Do not reproduce



What to collect? (I)

Vesicles:

- Use scalpel or 26 gauge needle to unroof vesicle – skin of roof goes to collection tube
 - Scrape base of vesicle with blunt edge scalpel, or wooden applicator and apply to microscope slide
- Lightly apply EM grid, shiny side down, against lesion. Repeat (X2) using more or less pressure.
- Repetitively touch a microscope slide to the lesion (touch prep)
- Allow slide, and grids to air dry for 10 minutes. Store in slide holder, and grid box, respectively

What to collect? (II)

Biopsy: 3.5-4 mm punch biopsy, bisect lesion, or obtain 2 biopsies.

1 sample in specimen collector tube (keep cold/frozen)


1 sample in formalin (not frozen).

Real-Time PCR assay: E9L-Vaccinia detection (Non-variola Eurasian orthopoxvirus assay)

- Samples are tested using primers and probe designed to detect Eurasian Orthopoxvirus other than variola
 - Potential human diseases detected:
 - **Vaccinia** **
 - Cowpox (Zoonotic disease of European origin)
 - Monkeypox (Zoonotic disease of central Africa)

Draft 1/20/03 Do not reproduce

Example of calibrated, R-T PCR data



Draft 1/20/03 Do not reproduce

Sensitivity of E9L vaccinia R-T PCR assay during validation at LRN labs

- **16/16 labs detected equivalent of 100 pfu vaccinia from dried, touch-prep slide...very sensitive!**

Draft 1/20/03 Do not reproduce

**If smallpox were to re-emerge,
test for vaccinia AE's would be
modified to become test for
variola virus DNA**

- **Alternate primer supplied...real-time PCR test otherwise essentially the same**

Vaccinia/orthopox and variola: real time PCR assays

- E9L: VAC, MPX, CPV(**TET**); variola (**FAM**) and ABI 7700
 - Essential gene of poxviruses
 - 16S control for inhibitors
 - Can be used to detect vaccinia (adverse event monitoring); use **TET** portion
 - Can be used to detect variola use **FAM** and **TET** portions (100% specificity in blind evaluation of 451 specimens)

Draft 1/20/03 Do not reproduce

Orthopoxvirus generic R-T PCR test

- Different E9L gene primer/probe targets *all* human pathogenic orthopoxviruses
- Currently being validated in LRN test sites as an “environmental sample” test
- Could have potential role in identification of non-vaccinia orthopoxvirus infections
- Sensitivity expected to be as high as test for detection of vaccinia

Draft 1/20/03 Do not reproduce



Draft 1/20/03 Do not reproduce




Draft 1/20/03 Do not reproduce




Draft 1/20/03 Do not reproduce




Draft 1/20/03 Do not reproduce



Draft 1/20/03 Do not reproduce





Draft 1/20/03 Do not reproduce



Draft 1/20/03 Do not reproduce

Negative Stain Electron Microscopy




vaccinia

**~1/2 hour per sample (for
experienced microscopist)**


Orthopoxvirus cell culture isolation

HEK 293




HeLa




Uninfected, 24 hpi



BSC-40



Variola (Ethiopia 17 R14-1X-72) infected, 24hpi



Draft 1/20/03 Do not reproduce

Selected issues relevant to implementation of pre-event vaccinia diagnostics

- Transition research/reference laboratory methods to routine patient care
 - Distribute test protocols
 - Test validation in multiple high-tech formats
 - Evaluation/validation of tests under 'field' conditions
 - Vaccinate testing personnel
 - QA/QC, Proficiency testing
- Reagents production/dispersal
 - Approval to use
- Development of vaccinia test algorithms (in conjunction with clinical algorithms)
- Involving private sector in rapid test development, licensure, deployment

Where to send orthopox specimens?

- Suspect **vaccinia adverse events** specimens that require identification of vaccinia go to closest **Laboratory Response Network (LRN)** laboratory.
 - Contact your State Public Health Lab Director for shipping address
 - All state and regional LRN labs can do real-time PCR for identification of vaccinia in AE's
- Specimens from persons with high suspicion of smallpox dx: Refer to *Rash, Vesicular Disease Algorithm*. Specimens go both to selected LRN with smallpox surge potential (contact CDC) and CDC simultaneously.

Draft 1/20/03 Do not reproduce



Specimen transport

- Standard dx specimen shipping guidelines available (subject to change):
<http://www.bt.cdc.gov/labissues/PackagingInfo.pdf>
- Serum, if collected, should be refrigerated and shipped
 - If spun and separated on site, freeze
- Formalin fixed material should be shipped at room temperature DO NOT FREEZE
- EM grids should be shipped at room temperature
- All other virus containing material should be stored and shipped frozen, unless it will be overnight shipped, then room temperature or refrigerated
- Keep all virus containing material out of direct sunlight

Past and Future

- What about the Past (when low tech worked)?
 - During smallpox epidemics *clinical* dx drove immediate medical response
 - Electron microscopy
 - Gel-diffusion serologic antigen detection
- Future? (dx development evolving rapidly)
 - Additional PCR targets (e.g., HA gene)
 - Antigen capture and DFA
- Post-event, expectations for dx would change

For More Information

- CDC Smallpox website
www.cdc.gov/smallpox
- National Immunization Program website
www.cdc.gov/nip