






Agents of Bioterrorism





Objectives

- Apply appropriate biosafety practices when working with potentially dangerous organisms
- Recognize agents that could be used for bioterrorism and perform selected tests to rule them out



Objectives

- Refer potentially dangerous specimens or organisms to an appropriate lab for confirmation
- Develop an Emergency Preparedness Plan for your laboratory.



Level A Laboratory: Definition

- **BSL₂-2 Laboratory with a certified Class II biological safety cabinet**
 - BSL-1 microbiology practices plus
 - Directed by competent scientists
 - Personnel specifically trained in handling pathogenic agents



Level A Laboratory: Definition

- BSL-2 Laboratory with a certified Class II biological safety cabinet
 - Physical containment practices to minimize infectious aerosols
 - PPE



Role of the Level A Laboratory

- Rule out critical biological agents
- Refer to higher level laboratory



Bioterrorism Agents: Laboratory Risk

<u>Agent</u>	<u>BSL</u>	<u>Laboratory Risk</u>
B. anthracis	2	low
Y. pestis	2	medium
Brucella spp.	2/3	high
F. tularensis	2/3	high
Botulinum toxin	2	medium
Smallpox	4	high
Viral Hemorrhagic fever	4	high



Francisella tularensis

Tularemia



Francisella tularensis

“A Rose by Any Other Name”

- Plague-like disease in rodents (California)
- Deer-fly fever (Utah)
- Glandular tick fever (Idaho and Montana)
- Market men's disease (Washington, DC)
- Rabbit fever (Central States)
- O'Hara's disease (Japan)





Tularemia

- Contagious --- no
- Infective dose --- 10-50 organisms
- Incubation period --- 1-21 days (average=3-5 days)
- Duration of illness --- ~2 weeks
- Mortality --- treated : low
untreated: moderate
- Persistence of organism ---months in moist soil
- Vaccine efficacy --- good, ~80%



Level A Procedures


Francisella tularensis

- This is a dangerous, highly virulent organism and it should not be manipulated at the bench.
- Gram stain
- Growth characteristics in broth
- Growth characteristics in agar



Francisella tularensis

- **Gram stain**
 - Poorly staining, tiny Gram-negative coccobacilli





Francisella tularensis

Growth Characteristics

- Fastidious, requires cysteine for robust growth: Cysteine Heart Agar (CHA) is ideal
 - Enriched chocolate agar : 9% sheep blood + cysteine
 - Not part of Level A routine procedures



Francisella tularensis **Growth Characteristics**





Gram Negative Coccobacilli

■ Most likely

- Acinetobacter
- Actinobacillus
- H. aphrophilus
- Bordetella spp.
- Pasturella spp.

■ Least likely

- DF-3
- Brucella spp.
- Francisella spp.



Francisella tularensis

Technical Hints

If you see:

- Tiny, gram-negative coccobacilli from blood, lymph node aspirate, or respiratory specimens
- Blood isolates that grow slowly on chocolate agar but poorly on blood agar
- Robust growth in BCYE; requires cysteine

Refer



Yersinia pestis

Plague



Plague Epidemiology

- U.S. averages 13 cases/yr (10 in 1998)
- 30% of cases are in Native Americans in the Southwest. 15% case fatality rate
- Most cases occur in summer



Plague Epidemiology

- **Three Clinical Types:**
 - bubonic (infected lymph nodes)
 - septicemic (blood-borne organisms)
 - pneumonic (transmissible by aerosol; deadliest)



Yersinia pestis

Specimen Selection

- Specimen selection is important
 - Bubonic - **bubo** - **lymph node aspirate**
 - Septicemic - **blood** - **Obtain three sets 10-30 minutes apart**
 - Pneumonic
 - **Sputum/throat**
 - **Bronchial washings**



Yersinia pestis

Specimen inoculation

- **Inoculate routine plating media and make thin smear for DFA**
 - **Use Wayson only if DFA is unavailable**



Level A Procedures

Yersinia pestis

- Gram stain
- Wayson stain
- Growth characteristics on agar
- Growth characteristics in broth



Yersinia pestis **Gram stain**

- Small, gram-negative bipolar-stained coccobacilli




Must confirm by DFA and mouse inoculation




Yersinia pestis **Wayson Stain**

- Pink-blue cells with a closed safety pin look





***Yersinia pestis* in BHI Broth**





Yersinia pestis


Technical Hints

- Small gram-negative, poorly staining rods from blood, lymph node aspirate, or respiratory specimens
- Safety pin appearance in Gram, Wright, Giemsa, or Wayson stain

Refer



Yersinia pestis blood plate





Variola virus

Smallpox








Variola Smallpox virus

- Large DNA virus
- Dumbbell-shaped core
- Complex membranes

Refer





Level A Procedures Smallpox virus

- **Rule out chickenpox (PCR)!**
- **Specimen of choice is lesion material from pustules**
- **Contact your State Public Health Laboratory for guidance**



Hemorrhagic Fever Viruses



Ebola



Marburg



Hemorrhagic Fever Viruses

- Families Responsible for VHF:
 - Arenaviridae
 - Bunyaviridae
 - Filoviridae
 - Flaviviridae



Hemorrhagic Fever Viruses

■ Arenaviruses

- Argentine Hemorrhagic Fever
- Bolivian Hemorrhagic Fever
- Sabia Associated Hemorrhagic Fever
- Lassa Fever



Hemorrhagic Fever Viruses

- **Bunyaviruses**
 - Crimean-Congo Hemorrhagic Fever
 - Rift Valley Fever
 - Hantavirus Pulmonary Syndrome
 - Hemorrhagic Fever



Hemorrhagic Fever Viruses

■ Filoviruses

- Ebola Hemorrhagic Fever
- Marburg Hemorrhagic Fever



Hemorrhagic Fever Viruses

- **Flaviviruses**
 - Tick-borne Encephalitis
 - Kyasanur Forest Disease
 - Omsk Hemorrhagic Fever



Viral Hemorrhagic Fevers

- Contagious --- Moderate
- Infective dose --- 1-10 particles
- Incubation period --- 4-21 days
- Duration of illness --- 7-16 days
- Mortality ---variable
- Persistence of organism --- unstable
- Non-endemic in U.S.
- No vaccine



VHF Specimens

- Diagnosis is clinical, not laboratory
- No specimen accepted without prior consultation

Refer



Handling VHF Specimens

- Sample for serology - 10-12 ml
 - ship on dry ice
- Tissue for immunohistochemistry
 - formalin-fixed or paraffin block
 - ship at room temperature
- Tissue for PCR/virus isolation
 - ante-mortem, post-mortem; ship on dry ice
- ship serum cold or on dry ice in a plastic tube



Brucella spp.

Brucellosis



BRUCELLOSIS

- Zoonotic disease caused by any of 4 *Brucella* spp.: *abortus*, *melitensis*, *suis*, and *canis*
- Systemic infection characterized by an undulant fever pattern
- Relatively rare in the U.S. with approximately 100 cases/year



Brucellosis

(by year, United States, 1966-1996)

BRUCE LLOSIS — by year, United States, 1966-1996



After peaking at more than 300 cases in 1975, the number of brucellosis cases has declined and, for the last 10 years, has remained relatively stable at approximately 100 cases per year.



BRUCELLOSIS: HISTORY

- 1887 Bruce - Malta fever, *M. melitensis*
- 1897 Bang - cattle abortion, *B. abortus*
- 1914 Traum - sow, *B. suis*
- 1920 Evans, Meyer, Shaw - Brucella
- 1954 *B. suis*, first weaponized U.S. agent
- 1968 Carmichael - Beagles, *B. canis*



BRUCELLOSIS: TRANSMISSION

- **Ingestion**
 - The most common mode of transmission
- **Direct skin contact/puncture**
 - Occupational hazard for farmers, butchers and veterinarians
- **Aerosols**
 - Highly infectious



BRUCELLSIS

- Infective dose = 10 -100 organisms
- Incubation period = 5 days - > 6 months
- Duration of illness = weeks to months
- Fever, profuse sweating, malaise, headache and muscle/back pain
- No person to person transmission
- Mortality: < 5%
- Stable organisms



Brucella spp. Specimen Selection

- Serum
- Blood or bone marrow
- Tissue (spleen, liver)



Level A Laboratory Tests

Brucella spp.

- **Colonial morphology on SBA**
- **Gram stain morphology**
- **Oxidase**
- **Urea hydrolysis**



B. abortus, gram stain (x3200)





***Brucella* spp.**

Key Level A Lab Tests

- **Colonial morphology on SBA**
 - Fastidious
 - Visible growth may take 48 - 72 hrs
 - Small (0.5-1.0mm), convex, glistening
 - Non-hemolytic and non-pigmented



Brucella spp.

Key Level A Lab Tests

- **Oxidase-positive**
 - *B. melitensis* (100%)
 - *B. abortus* (96%)
 - *B. suis* (95%)
 - *B. canis* (72%)
- **Urea hydrolysis-positive**
 - *B. suis* & *B. canis* ~15 min
 - *B. abortus* & *B. melitensis* ~24hr



Brucella spp.

Other Ox⁺/Ur⁺ GN Oxidizers

- **Achromobacter grp B**
- **Acidovorax spp**
- **Agrobacterium spp**
- **EO-2/EO-3**
- **Flavobacterium spp**
- **Methylobacterium spp**
- **Ochrobactrum**
- **Pseudomonas spp**
- **Riemerella**
- **Roseomonas spp**
- **O-2**
- **Il-i**



Brucella spp.

Technical Hints

- Misidentified as *Moraxella sp.*
 - Clin Inf Dis 1993; 17:1068-9
- Reported as gram-positive cocci, mistaken for “slow-growing” *Staphylococcus sp.*



Brucella spp.

Review of Key Tests

- Tiny, faintly staining, gram-negative coccobacilli from blood or bone marrow
- Oxidase +
- Urease +

REFER



Clostridium botulinum

Botulism




BOTULISM

- Diagnosis of botulism is made clinically
- Health care providers suspecting botulism should contact their State Health Department

REFER



Laboratory Capacity for Botulinum Toxin Testing





FOODBORNE BOTULISM

- Infective dose: 0.001 $\mu\text{g/kg}$
- Incubation period: 18 - 36 hours
- Dry mouth, double vision, droopy eyelids, dilated pupils
- Progressive descending bilateral muscle weakness & paralysis
- Respiratory failure and death
- Mortality 5-10%, up to 25%



FOODBORNE BOTULISM

- Among 309 persons with clinically diagnosed botulism reported to CDC from 1975 to 1988:
 - Stool cultures for *C. botulinum*: 51% +
 - Serum botulinum toxin testing: 37% +
 - Stool botulinum toxin testing: 23% +
- Overall, at least one of the above tests was positive for 65% of all patients



Level A Procedures for Botulism Event

- Properly collected specimens are to be referred to designated testing laboratories
- Prior to the shipment of any botulism-associated specimen, the designated laboratory must be notified and approved by the State Health Department



Level A Procedures for Botulism Event

- Clinical specimens to be collected:
 1. Serum
 2. Gastric contents or vomitus
 3. Feces or return from sterile water enema
 4. Wound tissue



Botulism Biosafety Alert

- Botulism toxins are extremely poisonous
- Minute quantities acquired by ingestion, inhalation, or by absorption can cause death
- All materials suspected of containing toxin must be handled with CAUTION!

REFER



Bacillus anthracis

Anthrax



ANTHRAX

■ Three forms of human anthrax occur:

1. Cutaneous
2. Gastrointestinal
3. Inhalational



Cutaneous anthrax



Vesicle development, day 2



Eschar formation, day 4



Anthrax Lesion on Neck





Inhalational Anthrax

- Infective dose = 8,000 - 15,000 spores
- Incubation period = 1-6 days
- Duration of illness = 3-5 days
- Fever, malaise, and fatigue
- Short period of improvement = up to 2 days
- Abrupt respiratory distress...death <24hrs
- No person to person transmission




Inhalational Anthrax

- Infection of tissue
- Late in the course of infection
- Chest X-ray



Mediastinal LN, microcolonies of *B anthracis*, Giemsa stain





Anthrax: Specimen Selection

- **Inhalation: Sputum and Blood**
- **Cutaneous: Vesicles and Eschar**
- **Gastrointestinal: Stool and Blood**




Bacillus anthracis

Key Level A Lab Tests

- Gram stain
- Growth characteristics on agar
- Sporulation, in air
- Lack of motility
- Penicillin inhibition zone
- Capsule formation



B. anthracis **Colony on SBA**






75

***B. anthracis* on blood (left) and bicarbonate agar (right) plate
(demonstrating capsule formation)**



“STICKY” consistency of *B. anthracis*’ colony on SBA

■





***B. anthracis*, Gram stain demonstrating spores**





Bacillus anthracis

Presumptive Identification

- Gram-positive rod
- catalase-positive
- Aerobic spore formation



Bacillus anthracis

Presumptive Identification

- Spores are oval
- non-swelling of vegetative cell
- ground glass colony morphology
 - *B. anthracis* (non-motile)
 - *B. cereus*
 - *B. cereus* var *mycoides* (non-motile)
 - *B. thuringiensis*



Bacillus anthracis

Presumptive Identification

- Nonmotile: *B anthracis* and *B cereus* var *mycoides* (and *B. megaterium*)
- Nonhemolytic, penicillin inhibition zone (15-20 mm)
- Capsule formation

REFER



CDC Laboratory Pearl

- The most common *Bacillus sp.* submitted to CDC to r/o *B. anthracis* are non-motile *B. megaterium*
- A non-motile, *Bacillus sp.*, is recovered from a blood culture. What is one of the fastest procedures to r/o anthrax ?



Acknowledgements

- Dr. Robbin Weyant
- Dr. Susan Maslanka
- Dr. May Chu
- Dr. Tom Ksiazek
- Dr. Peter Turnbull
- Dr. Edward Ewing